

## DATA SHEET

## Water purification system for autoclaves and labware washers Compact Sterifeed

Parameters	Sterifeed
Catalogue No.	CT-1002
TDS rejection rate	≥97%
Particles >0.22 μm	<1/ per mL
Nominal flow, pure water (to storage tank)	10 L/h
Deionization module	-
Display	Monochrome LCD
TOC monitor	-
Temperature range	15–30 °C
Relative humidity	20–80%
Mains voltage	200–240 VAC
Power	130 W
Fuse	230 V, 3.15 A
Width*Depth*Height	23*30*41 cm
System weight	12 kg
Operating weight	15 kg
Tank	Pressurized water storage tank required, not included

## **Feedwater properties**

Tap water
0.4–4.0 bar
<1500 µS/cm
5–35°C
4–10
<5
<0.1 ppm as CaCO <sub>3</sub>
<0.05 ppm as CaCO <sub>3</sub>
<0.05 ppm as CaCO <sub>3</sub>
<1 ppm
<+0.2
<2000 ppb



ADRONA SIA Kalnciema st. 209 Riga, LV-1046, LATVIA Phone: +371 67551894 Fax: +371 67551976 info@adrona.lv www.adrona.eu